



# Biopor<sup>®</sup>

## *Porous Polyethylene Implant*


A network of open, interconnecting pores supports fibrovascular ingrowth and integration with surrounding tissue. Strong, yet flexible, implants are easily contoured.

A unique water-soluble coating makes Biopor<sup>®</sup> implants easier to insert into tissue pockets. The coating reduces “tissue grab” and minimizes debris accumulation within the pores.

## COATED (A WATER-SOLUBLE AOC® COATING)


### Biopor® Sheets - AOC® Coated

	A	B	C
CSH-0554	0.45 mm	50 mm	40 mm
CSH-0558	0.45 mm	50 mm	80 mm
CSH-0854	0.8 mm	50 mm	40 mm
CSH-0858	0.8 mm	50 mm	80 mm
CSH-1554	1.5 mm	50 mm	40 mm
CSH-1558	1.5 mm	50 mm	80 mm
CSH-3054	3 mm	50 mm	40 mm
CSH-3058	3 mm	50 mm	80 mm




### Orbital Spheres - AOC® Coated

	Diameter
CSP-0016	16 mm
CSP-0018	18 mm
CSP-0019	19 mm
CSP-0020	20 mm
CSP-0021	21 mm
CSP-0022	22 mm




### Round Chin Implants, 2-piece, titanium connector - AOC® Coated

		A	B	C
CCR-4064	Small	60 mm	38 mm	4 mm
CCR-6064	Medium	60 mm	40 mm	6 mm
CCR-9064	Large	60 mm	43 mm	9 mm




### Square Chin Implants, 2-piece, titanium connector - AOC® Coated

		A	B	C
CCS-4064	Small	60 mm	38 mm	4 mm
CCS-6064	Medium	60 mm	40 mm	6 mm
CCS-9064	Large	60 mm	43 mm	9 mm




### Malar Implants - AOC® Coated

		A	B	C
CML-3052	Small Left	50 mm	20 mm	3 mm
CMR-3052	Small Right	50 mm	20 mm	3 mm
CML-4552	Medium Left	50 mm	20 mm	4.5 mm
CMR-4552	Medium Right	50 mm	20 mm	4.5 mm
CML-6052	Large Left	50 mm	20 mm	6 mm
CMR-6052	Large Right	50 mm	20 mm	6 mm



### Malar Shell Implants - AOC® Coated


		A	B	C
CSL-3052	Small Left	45 mm	23 mm	3 mm
CSR-3052	Small Right	45 mm	23 mm	3 mm
CSL-4052	Medium Left	49 mm	25 mm	4 mm
CSR-4052	Medium Right	49 mm	25 mm	4 mm
CSL-5052	Large Left	52 mm	25 mm	5 mm
CSR-5052	Large Right	52 mm	25 mm	5 mm



## UNCOATED


### Biopor® Sheets - Uncoated

	A	B	C
USH-0554	0.45 mm	50 mm	40 mm
USH-0558	0.45 mm	50 mm	80 mm
USH-0854	0.8 mm	50 mm	40 mm
USH-0858	0.8 mm	50 mm	80 mm
USH-1554	1.5 mm	50 mm	40 mm
USH-1558	1.5 mm	50 mm	80 mm
USH-3054	3 mm	50 mm	40 mm
USH-3058	3 mm	50 mm	80 mm




### Orbital Spheres - Uncoated

	Diameter
USP-0016	16 mm
USP-0018	18 mm
USP-0019	19 mm
USP-0020	20 mm
USP-0021	21 mm
USP-0022	22 mm




### Round Chin Implants, 2-piece, titanium connector - Uncoated

		A	B	C
UCR-4064	Small	60 mm	38 mm	4 mm
UCR-6064	Medium	60 mm	40 mm	6 mm
UCR-9064	Large	60 mm	43 mm	9 mm




### Square Chin Implants, 2-piece, titanium connector - Uncoated

		A	B	C
UCS-4064	Small	60 mm	38 mm	4 mm
UCS-6064	Medium	60 mm	40 mm	6 mm
UCS-9064	Large	60 mm	43 mm	9 mm




### Malar Implants - Uncoated

		A	B	C
UML-3052	Small Left	50 mm	20 mm	3 mm
UMR-3052	Small Right	50 mm	20 mm	3 mm
UML-4552	Medium Left	50 mm	20 mm	4.5 mm
UMR-4552	Medium Right	50 mm	20 mm	4.5 mm
UML-6052	Large Left	50 mm	20 mm	6 mm
UMR-6052	Large Right	50 mm	20 mm	6 mm



### Malar Shell Implants - Uncoated

		A	B	C
USL-3052	Small Left	45 mm	23 mm	3 mm
USR-3052	Small Right	45 mm	23 mm	3 mm
USL-4052	Medium Left	49 mm	25 mm	4 mm
USR-4052	Medium Right	49 mm	25 mm	4 mm
USL-5052	Large Left	52 mm	25 mm	5 mm
USR-5052	Large Right	52 mm	25 mm	5 mm



ORDER [orders@ceremed.com](mailto:orders@ceremed.com)  
310-815-2125

